

## Blu-DNA Stain Buffer (Cat. : BS-001)

[\(View DNA on electrophoresis gels with natural light\)](#)



### Description:

Blu-DNA stain buffer (1,000x) was developed to stain DNA or RNA on electrophoresis gels, detectable in the visible light ranges. Stored at room temperature, stable for 1 year.

**Detection:** Visible light/ nature light.

**Sensitivity:** Detect down to 100ng of DNA per band under nature light.

**Toxicity:** Non-mutagenic, Non-hazardous.

**Compatible:** Suitable for downstream cloning and sub-cloning applications.

### Protocol-1 for staining DNA during electrophoresis:

1. Prepare 1x Blu-DNA stain buffer in 0.5% - 2% of agarose gel solution containing 1x TAE/TBE and 0.5%-2% (w/v) agarose in a glass flask. (e.g. add 100ul of Blu-DNA Staining Buffer into 100ml of 1x TAE-agarose solution)
2. Heat in the microware until the solution is completely clear and no small floating particles are visible (about 2-3 minutes),
3. Pour the gel solution into a gel tray. After the agarose gel has solidified. The gel become blue and you can perform electrophoresis.
4. Mix 10ul of DNA samples (>100ng) with 2ul of 6x DNA loading buffer containing 1x TE buffer and 30% glycerol. Load samples into electrophoresis gel wells and run gel as usual.
5. (Optional) Place the gel into a water container after the electrophoresis if the DNA bands are weak. Gentle agitations 3 times for 5 minutes with distill water.  
(The DNA Bands are best visualized when viewed against a white background; light blue or, even better, on a light box)

### Protocol-2 for staining DNA after electrophoresis:

1. Add 100ul Blu-DNA stain buffer into 100ml 1x TAE or TBE running buffer to obtain 1x Blu-DNA stain buffer.
2. After electrophoresis, place gels in 1x Blu-DNA stain Buffer and gentle agitation for 15-30 minutes at room temperature or overnight at 4°C until the blue bands are visible on gel.
3. (Optional) Place the gel into a water container if the gel is over-stained. Gentle agitations 3 times for 5 minutes with distill water.  
(The DNA Bands are best visualized when viewed against a white background; light blue or, even better, on a light box)

**Benefits:** This method primarily eliminates the damage of DNA by uv irradiation and by Ethidium Bromide mutation. DNA isolated from blue stained gels should transform frozen competent E. coli (XL1-Blue and DH5) cells and obtain 20-50 folds more efficiently than Ethidium Bromide isolated DNA.

Contact us, Phone: 514-702 7702 Fax: 514-254 5356 Web: [www.zmtechscience.com](http://www.zmtechscience.com)

ZmTech Scientific Inc. endeavors to assist clients based on the highest level of customer service, competitive pricing and customer satisfaction. Our mission is:

Convenience, Speed, Safety and Economy.