

**Product Information:****Complete genomic DNA removal kit (Cat: DR-001)**

Kit contains: Reagent-A (1ml), reagent-B (400ul); 10x digestion buffer (400ul);
enzyme inactivation buffer (200ul);

Catalog: DR-001

Stored at -20°C.

Size: 200 reactions

Research use only, not for diagnostic use.

Description:

Complete genomic DNA removal kit contains a combination of nuclei acid enzymes and detergents that rapidly hydrolyze the entire genomic DNA, ds/ss-DNA, primer-dimers and oligos to mononucleotides, suitable for preparing DNA-free solutions (RNA solution, Taq DNA polymerase, PCR mastermix), used for PCR, RT-PCR, real-time PCR, Northern, microarray, RNA microinjection, RNA library construction, in vitro RNA synthesis and protein purification;

Benefits and features:

- In-tube digestions, no loss for your samples, no column involved for cleanup.

Protocol for the removal of genomic DNA contamination in RNA preparations

1. Mix 1ul of the reagent-A and 1ul of 10x digestion Buffer to the 18ul of RNA solution (up to 300ng/ul).

Gently mix by pipetting up and down.

2. Incubate reaction at 37 °C for 15-30 minutes.
3. Add 1ul of enzyme inactivation buffer and incubate reaction at 65 °C for 10 minutes.
4. Add 2ul of reagent-B into mixture. Incubate for 5 minutes at room temperature. Gently mix by pipetting up and down
5. Centrifuge 13,000 x g for 5 minute to pellet.
6. Transfer the supernatant (DNA-free RNA) to a new tube and discard the pellet (salts and enzymes).

The RNA is now ready to use for RT-PCR and other analysis.

Protocol for the removal of genomic DNA contamination in Taq DNA polymerase or PCR mastermix

1. Mix 5ul of the reagent-A and 2ul of 10x digestion Buffer to the 50ul of 2x PCR mastermix or Taq DNA polymerase (up to 500U). Gently mix by pipetting up and down.
2. Incubate reaction at 37 °C for 20-30 minutes.
3. Inactivate reaction at 95 °C for 3 minutes.

The PCR mastermix or Taq DNA polymerase is now ready to use for PCR and other analysis.

Tech tips: Not use for DNase I footprinting or random fragment analysis since the entire genomic DNA is completely digested.

Precautions and Disclaimer:

This product and procedure described are intended for R&D use only. Purchase of this product does not convey a license to perform any patented process.

Contact us,

Phone: 514-702 7702 Fax: 514-254 5356 Web: www.zmtechscience.com Email: order@zmtechscience.com (For ordering)