

**E-Chemiluminescence Reagents (Cat. E208082)****Product Information:**

- Contents:** Reagent-A (50mL) and Reagent-B (50mL),
- Catalog Number:** E208082
- Size:** 100 ml (1000 cm²)
- Storage Conditions:** Stored at 2-8°C
- Description:** Design for the detection of antibodies conjugated to Horseradish Peroxidase (HRP) in western blotting, suitable for both PVDF and NC membranes. Reagent A contains the Luminol and ECL enhancer. Reagent B contains the Hydrogen Peroxide (H₂O₂) and buffer stabilizer. The working solution (A&B mixture) is stable up to 12 hours at ambient temperature.

For research use only.

Procedure:

1. Prepare 1x working solution: mix **500ul reagent-A and 500ul reagent-B** for a mini-gel membrane (8 cm x 10 cm). The final volume of detection reagent mixture is around 0.1 ml/cm².
2. Place the membranes (protein side up) on a clean surface. Drain off the excess wash buffer.
3. Pipette the 1x working solution onto the membrane. The solution should cover the entire surface of the membrane.
4. Incubate for 5 minutes at room temperature without agitation.
5. Chemiluminescent detection:
Drain off excess working solution and place the blots (protein side up) on a clean surface.
Directly expose the membrane on a chemiluminescent / fluorescent imager or wrap up the membrane for x-ray film development.

Tech Tips:

1. Working solution mixture is stable up to 12 hours at ambient temperature.
2. Optimization range of primary antibody: 1/3,000-1/5,000.
3. Optimization range of secondary antibody (HRP-labeled): 1/30,000-1/50,000.

Precautions and Disclaimer:

This product and procedure described are intended for R&D use only. Purchase of this product does not convey a license to perform any patented process.

Contact us,

Phone: 514-702 7702 Fax: 514-254 5356 Web: www.zmtechscience.com

Email: order@zmtechscience.com (For ordering)