ECL Reaction Enhancer Buffer (Cat. E215060)

Product Information:

Contents: 2x ECL reaction enhancer

Catalog Number: E215060 Use for research only

Size: 250 mL Stored at 2-8°C

Description:

ECL reaction enhancer buffer is ready-to-use 2x ECL diluent, design for effectively increasing the ECL sensitivity (more than 3-10 folds) and boosting the protein signal intensity (light emission) in the complex reaction between peroxidase (HRP/AP), oxidant (H_2O_2) and luminol, suitable for most x-ray film or imager detection systems.

Procedure:

- 1. Prepare membrane after the incubation of secondary (2nd) antibody.
- 2. Prepare the chemiluminescence substrates (ECL) and ECL reaction enhancer buffer
 - 1) Mix equal volumes of ECL reagent-A and reagent-B (1:1)
 - 2) Add the same volume of 2x ECL reaction enhancer into the ECL mixture. The final volume of 1x working solution is around 0.1 ml/cm².
- 3. Wash membrane twice with washing buffer (1x TBST/PBST) for 5 minutes.
- 4. Place the membranes (protein side up) on a clean surface. Drain off the excess wash buffer.
- 5. Pipette the 1x working solution onto the membrane. The solution should cover the entire surface of the membrane.
- 6. Incubate for 5 minutes at room temperature without agitation.
- 7. Chemiluminescent detection:

Drain off excess working solution and place the blots (protein side up) on a clean surface.

Directly expose on a chemiluminescent / fluorescent imager or wrap up the membrane for x-ray film development.

Precautions and Disclaimer: This product and procedure described are intended for R&D use only. Purchase of this product does not convey a license to perform any patented process.

Contact us,

Phone: 514-702 7702 Fax: 514-254 5356 Web: www.zmtechscience.com Email: order@zmtechscience.com (For ordering)

ZmTech Scientific Inc. endeavors to assist clients based on the highest level of customer service, competitive pricing and customer satisfaction. Our mission is: Convenience, Speed, Safety and Economy.