

**Fast Blood Genomic DNA Extraction Kit (Cat. G217080)****Product Information:**

Contents: Reagent-A (50mL), Reagent-B (25mL), Reagent-C (2.5mL), Reagent-D (25mL), DNA precipitation solution (50mL)

Catalog Number: G217080

Size: 100 extractions

Storage Condition: stored at -20°C

For research use only

Description:


High-quality total genomic DNA (gDNA) is routinely extracted from anti-coagulated whole blood (EDTA or heparin), saliva, buccal swabs, or soft tissues, which might contain the bacteria, fungus, virus (HIV, HBV). The whole process in this kit provides a quick, easy and reliable method designed to remove the red blood cells (RBC), lysate all the nucleated cells (leukocytes), micro bacterial, fungus and viruses. The high-quality extracted genomic DNA is ideal for most PCR, southern blots, Microarray, restriction digestion and other diagnosis analysis.

Key features:

1. Obtain the highest yield and integrity of total genomic DNA from various blood sources within 1 hour.
2. No filter column or vacuum filtration is required, able to avoid the loss of small genomic DNA during extraction procedure.
3. No Phenol-Chloroform extraction and without using harmful chemicals.
4. High quality extracted DNA without RNA contamination, suitable for most downstream applications.

Procedure:

1. Add 500ul **Reagent A** into a 1.5 ml centrifuge tube contained 500ul anti-coagulated whole blood (EDTA or heparin). (Or 1:1 mix for large volume of blood).
Simple vortex and incubate for 5 minutes at room temperature.
2. Centrifuge at 12,000xg for 5 min to pellet the whole cells.
3. Discard the supernatant and add 225ul **reagent-B** to resuspend the pellet.
4. Simple vortex and incubate for 10 minutes at 37°C.
5. Add 25ul **reagent-C** into the sample tube and continue incubate for 10 minutes at 60°C.
6. Place the tube in a PCR machine (or water bath or block) and incubate at 95°C for 5 minutes.
7. Centrifuge at 12,000 xg for 2 minutes at 4°C and transfer the lysate supernatant into a clean tube.
8. Add 250ul **reagent-D** into the lysate supernatant and mix by gently vortexing.
9. Centrifuge at 5,000 xg for 5 minutes at 4°C and transfer the lysate supernatant into a clean tube.
10. Add 500ul **DNA precipitation solution** (Cat. #: PS-01D) into the lysate supernatant.
11. Mix well by pipetting up and down several times and centrifuge at 12,000 xg for 10 minutes at 4°C.
12. Aspirate liquid and simply wash pellet with 500ul 80% ethanol for 2 times by filling the tube with the ethanol solution and then decanting to discard the solution (don't resuspend the DNA pellets).
Centrifuge at highest speed (∞ 12,000 xg) for 10 minute at 4°C if the pellets are resuspended.
13. Air-dry pellet for 5-10 minutes and dissolve DNA in 30ul TE buffer or distilled water.
14. Centrifuge at 12,000 xg for 10 minutes at 4°C prior to measure the DNA concentration with 260/280nm. Store the genomic DNA solution at -20°C.
15. Pipette 0.5-1.5ul DNA solution into a 20ul PCR Master Mixture and run PCR/Real-Time PCR at thermal cycler.

**Flow Chart of genomic DNA Extraction:****Precautions and Disclaimer:**

This product and procedure described are intended for R&D use only. Purchase of this product does not convey a license to perform any patented process.

Contact us, Phone: 514-702 7702 Fax: 514-254 5356, Web: www.zmtechscience.com Email: order@zmtechscience.com