

**Product Information:****Fast Lysis-PCR Genotyping Kit**

Kit contains: Reagent-A (10ml); Reagent-B (500ul); Reagent-C (10ml); DNA precipitation solution (10ml);
Catalog: GT-001D
Sizes: 100 extractions
Storage : -20°C

Description:

Fast lysis-PCR genotyping reagents contain a combination of enzyme(s), detergents, and PCR required reagents that lysate the mouse/rat tissues (tails, ears, yolk sacs) within 30 minutes and directly used to run PCR reactions without any further DNA purification.

Procedure:

1. Prepare 0.2-0.3cm mouse tail biopsy sample in a 0.5ml or 1.5ml microcentrifuge tube.
2. Add 100ul reagent-A and 5ul reagent-B into the sample tube. Incubate at 60°C for 10 minutes.
3. Place the tube in a PCR machine (or water bath or block) and incubate at 95°C for 20 minutes.
4. Centrifuge at 12,000 xg for 2 minutes at 4°C and transfer the lysate supernatant into a clean tube.
5. Add 100ul reagent-C into the lysate supernatant and mix by vortexing.
6. Pipette 2-10ul lysates into a 20ul PCR mastermix and run PCR/Real-Time PCR at thermal cyclers.


Suggested PCR Protocol:**I. Preparation of PCR Master Mix**

for a single reaction (total volume: 20uL) in a 0.2 or 0.5mL microtube.

| Component | Volume (μL) | Final Concentration |
|----------------------|-------------|---------------------|
| Lysate supernatant | 2-10 | determined by user |
| 2x PCR mastermix | 10 | 1x |
| Forward primer (5μM) | 1 | 200nM |
| Reverse primer (5μM) | 1 | 200nM |
| PCR grade water | up to 20 μL | |

II. Setup typical thermal cycling parameters

| | | |
|-------------------------|------|---|
| Enzyme activation step: | 95°C | 5 minutes |
| 25-40 cycles: | | |
| Denaturation | 95°C | 30 seconds |
| Annealing | X°C | 30 seconds dependent on T _m of primers |
| Extension | 72°C | 30 seconds (1min per kb amplicon) |
| Hold at 4-8°C | | |

III. Figure:



PCR analysis of genomic DNA extracts from mouse tail biopsy at 0.2-0.3 cm, using Fast Lysis-PCR genotyping kit (GT-001). Genomic DNA in 2ul lysate supernatant was directly amplified in a 20ul PCR reaction for 35 cycles.

Technical Tips:

1. Ensure the tissue samples (tail, ear or yolk sac) are completely submerged in the reagents.
(Use roughly proportional volume of reagents for different sized samples)

| | mouse tail (0.2-0.3cm) | ear(0.2cm) | yolk sac (E8.5) |
|-----------------|------------------------|------------|-----------------|
| Reagent-A (ul): | 100 | 80 | 50 |
| Reagent-B (ul): | 5 | 4 | 2.5 |
| Reagent-C (ul): | 100 | 80 | 50 |

2. Tissues will not completely digested at the end of the incubation. This is normal and not influence PCR performance. Undigested tissues were stored at -20°C, used for other applications or further genotypings. If necessary, longer incubation or overnight at step 2. can completely digest the tissue sample and no loss of efficacy.
3. Lysates are stable at 4°C for at least 6 months and at room temperature for 1-2 weeks. For longer-term storage at room temperature, it is necessary to further purify the crude lysates, using the following procedure: add NaCl to a final concentration of 250mM, and then add 0.7 volume of pure ethanol to precipitate DNA, spin down DNA at 4°C for 3 minutes, discard supernatant and dissolve DNA in 50ul TE buffer (10 mM Tris-Cl, pH 7.5; 1 mM EDTA).
4. The kit may be used for other genomic DNA extractions from various animal tissues, hair shafts, toes or saliva.
5. Additional information: (An new technical procedure for DNA precipitation)

Procedure:

1. Prepare 0.2-0.3cm mouse tail biopsy sample in a 0.5ml or 1.5ml microcentrifuge tube.
2. Add 100ul reagent-A and 5ul reagent-B into the sample tube. Incubate at 60°C for 10 minutes.
If necessary, stroke the tissues with pestle can completely digest the tissue sample to obtain higher DNA yields.
3. Place the tube in a PCR machine (or water bath or block) and incubate at 95°C for 20 minutes.
4. Centrifuge at 12,000 xg for 2 minutes at 4°C and transfer the lysate supernatant into a clean tube.
5. Add 100ul reagent-C into the lysate supernatant and mix by vortexing.
6. Centrifuge at 3,000 xg for 5 minutes at 4°C and transfer the lysate supernatant into a clean tube.
7. Add 100ul DNA precipitation solution (Cat.#: PS-01D) into the lysate supernatant. Mix well and centrifuge at 12,000 xg for 10 minutes at 4°C. Aspirate liquid and simply wash pellet with 80% ethanol for 2 times (don't resuspend the DNA pellets). Air-dry pellet for 5-10 minutes and dissolve DNA in 50ul TE buffer or distilled water. Measure the DNA Concentration with 260/280nm and store DNA at room temperature or at 4°C.
8. Pipette 1-2ul lysates into a 20ul PCR Master Mixture and run PCR/Real-Time PCR at thermal cyclers.

Precautions and Disclaimer:

This product and procedure described are intended for R&D use only. Purchase of this product does not convey a license to perform any patented process.