

**Product Information****Inova D-gel solution (4x)****Catalog Number:** DG-001**Stored at room temperature****Description:**

**Inova D-gel solution** provides a new gel metric for analyzing the sizes of protein, nucleic acid (DNA/RNA/oligo) contains acrylamide and agarose in liquid phase.

- Very clear background
- easy handle and gel making.

**Protocol-I. using LB-001 loading buffer:**

1. Prepare the hot 30 ml of 1xTAE, TBE or Borate Buffer (>70oC) without **Ethidium Bromide** in a glass flask.
2. Add 10 ml of the D-gel solution into the glass flask, gentle shaking to complete dissolve the D-gel solution and pour into a gel tray. After the agarose gel is solidified, then perform electrophoresis.
3. Add 1-2 µl of the 6X fluo- DNA/RNA loading buffer (LB-001/LB-201) to 5-10 µl DNA/RNA samples. Mix thoroughly.
4. Load DNA/RNA samples and run the gel using your standard protocol.
5. View DNA/RNA bands using a **blue/UV light transilluminator** during or after electrophoresis.
6. Images can be taken using a blue light transilluminator or a UV transilluminator.

**Protocol-II. using GS-001 gel staining solution:**

1. Prepare the hot 30 ml of 1xTAE, TBE or Borate Buffer (>70oC) without **Ethidium Bromide** in a glass flask.
2. Add 10 ml of the D-gel solution and 2 ul of 30,000x Fluo-DNA/RNA gel staining solution (GS-001) into the glass flask, gentle shaking to complete dissolve the D-gel solution and pour into a gel tray. After the agarose gel has solidified you can perform electrophoresis.
3. Load DNA/RNA samples and run the gel using your standard protocol.
4. View DNA/RNA bands using a **blue/UV light transilluminator** during or after electrophoresis.
5. Images can be taken using a blue light transilluminator or a UV transilluminator.

Precautions and Disclaimer: This product and procedure described are intended for R&D use only. Purchase of this product does not convey a license to perform any patented process.