



ZmTech® Fluo-DNA Loading Buffer

Product Information:

Contains: 6x Fluorescent DNA dye loading buffer (RNase-free and Dnase-free, sterilized).
Catalog: LB-001, LB-001G, LB-001B, LB-001R, LB-001GB, LB-001BX, LB-001GR, LB-001GX
Sizes: 1.0ml
Storage Conditions: Store at -20°C, protected from light, stable for 1 year.

Order information (product series):

- **Cat.: LB-001** contains 6x fluorescent DNA dye without loading/tracking dye
- **Cat.: LB-001G** contains 6x fluorescent DNA dye with loading/tracking dye (**Orange G** : ~50bp)
- **Cat.: LB-001B** contains 6x fluorescent DNA dye with loading/tracking dye (**Bromophenol blue**: ~400bp)
- **Cat.: LB-001R** contains 6x fluorescent DNA dye with loading/tracking dye (**Cresol red**: ~1kb)
- **Cat.: LB-001GB** contains 6x fluorescent DNA dye & tracking dyes (**orange G / Bromophenol blue**:~50- 400bp)
- **Cat.: LB-001BX** contains 6x fluorescent DNA dye & tracking dyes (**Bromophenol blue/Xylene cyanol**:400-4kb)
- **Cat.: LB-001GR** contains 6x fluorescent DNA dye with tracking dyes (**Orange G/Cresol red**: ~50bp-1kb)
- **Cat.: LB-001GX** contains 6x fluorescent DNA dye with tracking dyes (**Orange G /Xylene cyanol**: ~50bp-4kb)

Features and Benefits:

1. **Safe:** Non-toxic and non-genotoxicity for waste disposal, directly into the wastewater systems or clean up with water and 70% ethanol. No need to expose your skin and eyes on UV light.
2. **Sensitive:** Detect down to 60 ng of DNA per band with blue light and 1-5 ng of DNA per band with UV light.
3. **Convenient:** Visualize DNA bands with a blue light transilluminator or U.V transilluminator during/after electrophoresis in TAE/TBE buffer.
4. **Effective:** Eliminating ultraviolet and EtBr induced mutation or cleavage of DNA fragments.
5. **Compatible:** Fluorescent DNA dye can be completely removed from nucleic acids by alcohol precipitation or Qiagen QIAquick Gel Extraction, suitable for downstream cloning/sub-cloning applications.

DNA Staining Protocol:

1. Prepare 0.5% to 2% of agarose gel solution in 1xTAE, TBE or Borate Buffer without **Ethidium Bromide** in a glass flask.
2. Heat in the microwave until the solution is completely clear and no small floating particles are visible (about 2-3 minutes). Pour the gel solution into a gel tray after cool down. After the agarose gel has solidified you can perform electrophoresis.
3. Add 1-2 µl of the 6X fluo- DNA loading buffer (LB-001) to 5-10 µl DNA sample or DNA marker/ladder. Mix thoroughly.
4. Load samples and run the gel using your standard protocol.
5. View DNA bands using a **blue light transilluminator** during or after electrophoresis.
6. Images can be taken using a blue light transilluminator or a UV transilluminator.

Technical Tips:

1. Loading 60ng-200ng DNA each lane is ideal for viewing DNA bands in agarose gels. Too much would affect the DNA migration.
2. Mix well the DNA samples or ladders with 6x fluo-DNA loading buffer by pipetting up and down for several times.
3. For removing Fluo-DNA dye from DNA samples by simple ethanol precipitation: add NaCl to a final concentration of 250mM, and then add 0.7 volume of pure ethanol to precipitate DNA, incubate on ice or -20°C for 20 minutes and spin down DNA at 4°C for 10 minutes, discard supernatant and dissolve DNA in TE buffer (10 mM Tris-Cl, pH 7.5; 1 mM EDTA).
4. **Loading/tracking dye information:**

Colour dye	0.5-1.5% agarose	2.0-3.0% agarose	CAS number	Cat. No
Xylene cyanol	10'000-4000 bp	750-200 bp	2650-17-1	Sigma X4126
Cresol Red	2000-1000 bp	200-125 bp	62625-29-0	Sigma 114480
Bromophenol blue	500-400 bp	150-50 bp		Sigma B8026
Orange G	<100 bp	---		Sigma O3756

Precautions and Disclaimer: This product and procedure described are intended for R&D use only. Purchase of this product does not convey a license to perform any patented process. **Contact us,** Phone: 514-702 7702 Fax: 514-254 5356 Web: www.zmtechscience.com