



## RNA Precipitation Solution (Cat. PS-02R)

### Product Information:

**Catalog Number:** PS-02R  
**Sizes:** 10 mL (2x), nuclease-free, sterile  
**Storage :** 2-8°C

### Description:

- The **RNA Precipitation Solution (2x)** is designed for rapid purification and concentration of RNA during RNA extractions. This precipitation solution offers major advantages over the ethanol or isopropanol precipitation methods in that it **selectly** precipitates RNA and does not precipitate DNA, protein or carbohydrates. It is able to remove genomic DNA, Proteinase K and Nuclease during RNA preparation. Moreover, it helps to remove inhibitors of cDNA synthesis and maximizes the yields of total RNA, specifically extracted from the small/tiny tissues or LCM samples. The RNA precipitation with >90% efficiency is accomplished in 5-10 minutes on ice. The concentrated RNA is suitable for most downstream applications: real time PCR, reverse transcription PCR, northern blotting, RNase protection, primer extension, hybridization, cloning, microarray assays or other RNA assays.

### Procedure:


- Add an equal volume of RNA Precipitation Solution (2x) into the RNA fragment solutions (cell/tissue lysates). Mix well and incubate at -20°C for 20 minutes or overnight.
- Centrifuge at 12,500 x g for 15 minutes at 4°C.
- Carefully aspirate liquids and simply rinse pellets with 2 volumes of ice cold 80% ethanol for 2 times without resuspending the RNA pellets. The RNA pellets will **not** be visible if the concentration is less than 20ng/ul.
- Dissolve the RNA pellets in 20 µl of **nuclease-free** H<sub>2</sub>O or TE if the RNA pellet is visible. Otherwise, use 10ul of **nuclease-free** H<sub>2</sub>O or TE.
- Measure RNA concentration using a spectrometer and store RNA at -80°C.

Note: This precipitation solution may not remove all DNA contamination from total RNA isolations. If the RNA is to be used for RT-PCR, it should be DNase I treated.

### Suggested Protocol for DNase I treatment:

- Mix 20ul RNA solution ( $\geq 10\mu\text{g}$  RNA) with 20ul of 2X DNase I Reaction Buffer (20 mM Tris-HCl, 5 mM MgCl<sub>2</sub>, 1 mM CaCl<sub>2</sub>, pH 7.6)
- Add 2 units of DNase I, mix thoroughly and incubate at 37°C for 10 minutes.
- Add 1 µl of 0.5 M EDTA (to a final concentration of 5 mM). Heat inactivate at 75°C for 10 minutes

### Figures:



The analysis of crude RNA extracted from MCF-7 cells in 1.2% Formaldehyde RNA gels:

- (A) Before using the RNA precipitation solution (PS-02R) showed gDNA (top band), 28S, 18S and 5S RNA bands (lower bands)  
(B) After using the RNA precipitation solution (PS-02R) only showed 28S, 18S and 5S RNA bands.

### Precautions and Disclaimer:

This product and procedure described are intended for R&D use only. Purchase of this product does not convey a license to perform any patented process.

Contact us, Phone: 514-702 7702 Fax: 514-254 5356 Email: [order@zmtechscience.com](mailto:order@zmtechscience.com) Web: [www.zmtechscience.com](http://www.zmtechscience.com)