



Western Blot Stripping Buffer (Cat. S238070)

Product Information:

Contents: Western blot stripping buffer
Catalog Number: S238070
Size: 500 mL
Storage Conditions: Store at room temperature

Description:

Western Blot Stripping Buffer (S238070) is suitable to remove the 1° antibodies, the secondary antibodies and various chemiluminescent substrates on PVDF membranes. No loss of target proteins, enable to detect multiple phosphorylated and nonphosphorylated proteins on the same membrane.

Features and Benefits:

1. No loss of transferred proteins; No odor; No toxic chemicals involving;
2. **Working at room temperature;**
3. More than 10-20 times of re-probing processes on the same membrane (PVDF).

Procedure:

1. Wash membrane twice for 5 minutes with washing buffer (1x TBST/PBST) to remove chemiluminescent substrates on membrane surface.
2. Incubate and agitate the membrane in 1x western blot stripping buffer (Cat.S238070) for 15-30 minutes at room temperature with shaking.
3. Wash membrane 3 times for 10 minutes with washing buffer (1x TBST/PBST).
4. Block membrane in Zmtech protein block-on (Cat.B208060) or other blocking buffer for 1 hour at room temperature.
5. Re-probed with first antibodies followed by goat anti-mouse/rabbit HRP/AP and detected by chemiluminescence substrates.

Precautions and Disclaimer:

This product and procedure described are intended for R&D use only. Purchase of this product does not convey a license to perform any patented process.

Contact us, Phone: 514-702 7702 Web: www.zmtechscience.com Email: order@zmtechscience.com (For ordering)

ZmTech Scientific Inc. endeavors to assist clients based on the highest level of customer service, competitive pricing and customer satisfaction. Our mission is: Convenience, Speed, Safety and Economy.