Simple Transfection Reagent (Cat: T3080)

Sizes: 500uL (5x 100ul) Research use only, not for diagnostic use.

Storage: -20°C or -80°C. • After thawing the solution can be stored at 4 °C for up to 2 months.

Description:

Simple transfection reagent is a ready-to-use polymer-based transfection agent, with high transfection efficiency and lowest toxicity to deliver DNA or siRNA into various primary/ mammalian cell lines.

Workflow Timeline:

Day 1:Seed cellsDay 2:(am/pm)Transfect Cells (am), After 2h-4h transfection - add the complete media (pm).Day 3:(pm) or moreObserve fluorescence, or by FACS analysis.

General Protocol:

Day 1: Cell Seeding: Split cells one day before transfection.

a. <u>6 well plate:</u> $2x10^5$ cells/well

b. <u>10 cm dish:</u> 1×10^6 cells/dish

The cell density should reach to the 40-90% confluency at the time of transfection.

Day 2: (Am) Prior to transfection, bring Simple-Transfection reagent (T3080) to room temperature and warm the serum-free media (DMEM w/o phenol red) or PBS to 37°C.

1. Preparation the complex: <u>Transfection agent (ul) + DNA(ug) (3:1)</u>.

* Dilute total plasmid DNA (ug)/siRNA in serum-free media (DMEM w/o phenol red) or PBS.

- * Add Simple Transfection agent to the diluted DNA/RNA. (3:1 ratio of T3050 (ul): DNA (ug))
- a. For <u>6 well plate:</u> 300ul serum-free media (DMEM)+ **3 ug of DNA**, vortex briefly to mix well, add **9 ul of Simple- Transfection reagent;**
- **b.** For <u>10cm dish:</u> 900uL serum-free media (DMEM)+ 9 ug of DNA, vortex briefly to mix well, add 27ul of Simple-Transfection reagent;
- 2. Immediately pulse vortex for 10 seconds. Spin down briefly.

3. Incubate 15 minutes at Room Temperature (RT) for forming the TransFection/DNA or

RNA complexes.(The formation of TransFection/DNA might make solution slightly cloudy).

4. Aspirate the cell culture medium from each well/dish and add the Transfection/DNA complexes to cells carefully drop-wise. Gently rock the plate to ensure even distribution. Do not swirl plate and not to dislodge the cells.

5. Return the well/dish to a 37 °C incubator with 5% CO₂.

Day 2: (Pm) Incubate at 37°C for 2h–4h, add appropriate volumes of complete cell culture medium to the well/ dish and continue to incubate the cells for 12-24 h. a. For <u>6 well plate:</u> add 2 mL complete medium/each well.

b. For 10 cm dish: add 9 mL complete medium/dish.

Day 3: (Pm) or more: Observe fluorescence, harvest cells, or perform your experiment by FACS analysis.

Precautions and Disclaimer: This product and procedure described are intended for R&D use only. Purchase of this product does not convey a license to perform any patented process.