



ZmTech Scientific Innovation and Development of Laboratory Accessories and Bio-Technologies

® • **ZmTech Mycoplasma PCR/RT-PCR Detection Kits (Cat. M208001/M208002)**

Kit contains: (All-IN-ONE PCR Mastermix)

1. Ready-to-use SYBR Green PCR Mixtures, including optimization of 1x PCR Buffer, Primer Sets, Taq DNA Polymerase, dNTP (dATP, dCTP, dGTP, dTTP), PCR Buffer Stabilizer, RNase/DNase inhibitors and SYBR Green I dye (Blue Cap)
2. Positive control (Yellow Cap) Stored at -20°C

Supplemental Requirements:

1. PCR/Real-time PCR devices
2. PCR reaction tubes or Glass Capillaries

Features and Benefits:

1. Easy: Samples are prepared directly from culture cell supernatant, frozen cells or cell lysates.
2. Safest: Samples directly get from frozen cells, not need to re-culture cells in incubator and handle in hood.
3. Sensitive: able to detect 1fg Mycoplasma DNA, equivalent to one or two genome copies of the 16S rRNA coding region.
4. Broadest: able to detect up to 117 species of Mycoplasma. (see P.S.*)

Protocol: (Sample can be prepared from culturing cells or frozen cells or cell Lysates)

1. Sample preparation A: 100ul supernatant from cell culture medium and heat at 95°C for 5 minutes.
 Sample preparation B: 1-10 x 10³ cells in 100ul PBS from frozen cells and heat at 95°C for 5 minutes. (Note*)
 Sample preparation C: 2ul cell Lysates in 100ul PBS and heat at 95°C for 5 minutes.
 (*Briefly centrifuge the heated samples for 5-10 seconds to pellet cellular debris before adding to the PCR mixture.)
2. Run reactions in the spectrofluorometric thermal cyclers (Suggested Protocols: *SP)

***SP 1.**

Run reactions in the Normal PCR thermal cyclers as following cycling parameters in a 25ul plastic tube.

Add 2.5ul Samples /2.5ul positive control / 2.5ul water (Negative)

into 22.5ul Ready-to-use Mixtures (Blue cap)

Total: 25ul running reaction into a PCR Tube

	Step:	Temperature:	Time:
1.	Denaturation (1 cycle)	94°C	10 minutes(Essential)
2.	Cycles (39 cycles):		
	Denature	94°C	30 seconds
	Annealing	55°C	30 seconds
	Extension	72°C	30 seconds
3.	Cool down to	4-8°C	End → Run 2% agarose gel at 120v for 30 minutes

***SP 2.**

Run reactions in the Real-Time PCR thermal cyclers as following cycling parameters in a 25ul plastic tube.

Add **2.5ul** Samples /2.5ul positive control / 2.5ul water (Negative)

into **22.5ul Ready-to-use Mixtures (Blue cap)**

Total: 25ul running reaction into a PCR Tube

Step:	Temperature:	Time:
1. Denaturation (1 cycle)	94°C	10 minutes (Essential)
2. Cycles (39 cycles):		
Denature	94°C	30 seconds
Annealing	55°C	30 seconds
Extension	72°C	30 seconds
3. Melting(Dissociation) Curve Analysis	(see instrument manufacturer guidelines)	

***SP 3.**

Run reactions in the Real-Time PCR thermal cyclers as following cycling parameters in a 20ul capillary tube.

Add **2 ul** Samples /2ul positive control / 2ul water (Negative)

into **18 ul Ready-to-use Mixtures (Blue cap)**

Total: 20ul running reaction into a capillary tube

Step:	Temperature:	Time:	Temperature Transition Rate
1. Denaturation (1 cycle)	94°C	10 minutes	20°C/second
2. Cycles (35 cycles):			
Denature	94°C	15 seconds (none)	20°C/second
Annealing	55°C	15 seconds (none)	20°C/second
Extension	72°C	8 seconds (none)	20°C/second
Detection	72°C	0 seconds (single)	20°C/second
3. Melting(1 cycle)	52°C	1 minute (none)	20°C/second
		(95°C, 0 Second, 0.15°C/second, con.)	
4. Cooling (1cycle)	40°C	30 seconds (none)	20°C/second

(Note*):

Simply centrifuge cells and re-suspend in 100ul PBS for removing DMSO in medium. After heated, the samples are vortex for 5 seconds and centrifuged at 13,000 rpm for 2 minutes at room temperature. This supernatant will be used as DNA template in PCR detection.

Contact us,

Phone: 514-702 7702 (Sam Lee) Fax: 514- 254 5356 (Montreal)

Email: ZmTech.services@hotmail.com or Sam.lee@zmtechsci.com (For ordering @ Technical supports)

Web: www.zmtechsci.com

Test Evaluation:

- A) A 167bp or 130bp DNA band was detected if running contaminated samples in 1.5-2% Agarose gel for 30 minutes at 120V.
- B) Approximately evaluate the density of Mycoplasma contaminations from representative RT-PCR amplification curves (Delta Rn vs Cycle)
- C) Determine some Mycoplasma species by melting(dissociation) peak Tm°C

Figure1. Gel Evaluation: a distinct 167bp or 130bp band indicates positive control or Mycoplasma contamination.




Figure2. Evaluate Mycoplasma genome copies in each PCR reaction and determine the density of Contamination.





Figure3.Using the dissociation-curve analysis (melting peak value) to determine Mycoplasma species in cell lines:



Species:	<i>M. Fermentans</i>	<i>M. .hyorhinis</i>	<i>M. orale</i>	<i>M. Salivarium</i>	<i>M. arginini</i>
Tm°C:	79.4°C (H)	79°C (H)	79.2°C(H)	77.6°C(H)	77.7°C(H)
Tm°C:	81.0°C (M)	81.0°C (M)	81.5°C(P)	79.5°C(M)	81.5°C(R)

P.S*: Mycoplasma Species Detected by

® •ZmTech Mycoplasma PCR/RT-PCR Detection Kit (Cat. M208001/M208002)

1	<u>Acholeplasma entomophilum</u>	41	<u>M. cricetuli</u>	81	<u>Mesoplasma entomophilum</u>
2	<u>Acholeplasma modicum</u>	42	<u>M. cynos</u>	82	<u>M.moatsii</u>
3	<i>A. laidlawii</i>	43	<i>M. dispar</i>	83	<u>M.mobile</u>
4	<u>A. oculi</u>	44	<i>U. diversum</i>	84	<u>M. molare</u>
5	<u>Anaerplasma anaerobium</u>	45	<u>Entomoplasma somnilux</u>	85	<u>M. mustelae</u>
6	<u>Asteroleplasma abastoclassicum</u>	46	<u>M. ellychniae</u>	86	<u>M. muris</u>
7	<u>Asteroleplasma varium</u>	47	<i>M. equirhinis</i>	87	<i>M. mycoides</i> subsp
8	<u>M. adleri</u>	48	<i>M. equigenitalium</i>	88	<i>M. neurolyticum</i>
9	<i>M. agalactiae</i>	49	<i>M.falconis</i>	89	<i>M.opalescens</i>
10	<i>M. agassizii</i>	50	<u>M. fastidiosum</u>	90	<i>M. orale</i>
11	<u>M. alkalescens</u>	51	<u>M. faucium</u>	91	<i>M. ovipneumoniae</i>
12	<u>M.alligatoris</u>	52	<u>M.felifaucium</u>	92	Ovine ureaplasmas
13	<u>M.anseris</u>	53	<u>M. felis</u>	93	<i>M.oxoniensis</i>
14	<i>M. anatis</i>	54	<i>M. fermentans</i>	94	<i>U.parvum</i>
15	<i>M. arginini</i>	55	<i>M. flocculare</i>	95	<i>M. penetrans</i>
16	<i>M. arthritis,</i>	56	<i>M.gallinarum</i>	96	<i>M.phovirninis</i>
17	<i>M. auris</i>	57	<i>M.gallinaceum</i>	97	<i>M. pirum</i>
18	<u>M. bovirhinis</u>	58	<u>M. gallisepticum</u>	98	<i>M. pneumoniae</i>
19	<u>M. bovirhinis</u>	59	<u>M.gallopavonis</u>	99	<i>M.primatum</i>
20	<u>Bovine Group 7</u>	60	<u>M. gateae</u>	100	<i>M. pulmonis,</i>
21	<i>M. bovis</i>	61	<i>M. genitalium</i>	101	<i>M. pullorum</i>
22	<u>M. bovoculi</u>	62	<i>M. glycyphilum</i>	102	<u>M. putrefaciens</u>
23	<u>M. buccale</u>	63	<i>M.gypis</i>	103	<i>M. salivarium</i>
24	<u>M. buteonis</u>	64	<i>M. hominis</i>	104	<i>M. SP1</i>
25	<u>M. californicum</u>	65	<i>M. hypopneumoniae</i>	105	<i>M. SP2</i>
26	<u>M. canadense</u>	66	<u>M. hyopharyngis</u>	106	<i>M. spermatophilum</i>
27	<u>M. canis</u>	67	<u>M. hyorhinis cultivar a</u>	107	<i>M. spumans</i>
28	<u>Capri</u>	68	<u>M. hyorhinis</u>	108	<i>Spiroplasma apis</i>
29	<u>Capripneumoniae</u>	69	<u>M. hyosynoviae</u>	109	<u>Spiroplasma citri</u>

30	<i>M. capricolum</i>)	70	M.iguanae	110	M. stumi
31	<i>M. Caviae</i>	71	M. imitans	111	<i>M. subdolum</i>
32	M. cavipharyngis	72	M. indiense	112	<i>M.sualvi</i>
33	M.citelli	73	<i>M. iners</i>	113	M. synoviae
34	M.cloacale	74	M. iowae	114	Ureaplasma urealyticum
35	M.collis	75	M.lagogenitalium	115	<i>U. canigenitalium</i>
36	M. columbinassale	76	<i>M. leopharyngis</i>	116	<i>M. verecundum</i>
37	M. columbinum	77	<i>M. lipofaciens</i>	117	M. yeatsii
38	M. columborale	78	M. lipophilum		
39	M. conjunctivae	79	M. maculosum		
40	M. corogypsi	80	M. meleagridis		

Other bacterial strains are not be Detected by

® **ZmTech Mycoplasma PCR Detection Kits (Cat. M208001/M208002)**

Actinomyces israelii	Chromatium vinosum	Haemophilus influenzae	Ornithobacterium rhinotracheale	Staphylococcus pneumoniae
Bacillus subtilis	Clostridium innocuum	Haemophilus parainfluenzae	Pasteurella avium	Staphylococcus pyogenes
Bacteroides nodosus	Clostridium ramosum	Helicobacter felis	Pasteurella gallinarum	Streptomyces ambofaciens
Borrelia burgdorferi	Corynebacterium variabilis	Lactobacillus cateniforme	Pasteurella multocida	Streptococcus mutans
Campylobacter hyointestinalis	E. coli	Lactobacillus fermentum	Propionibacterium acnes	Streptococcus pneumoniae
Candida albicans	Escherichia coli	Megasphaera elsdenii	Pseudomonas aeruginosa	Streptococcus pleomorphus
Chlamydia psittaci	Erysipelothrix rhusiopathiae	Mycobacterium tuberculosis	Salmonella Typhimurium	Thermus thermophilus
Chlamydia trachomatis	Eubacterium biforme	Neisseria gonorrhoeae	Staphylococcus aureus	
Chlamydia pneumoniae	Gardnerella vaginalis	Nocardia asteroides	Staphylococcus epidermidis	

Reference:

Agata Baczynska, Helle F Svenstrup, Jens Fedder, Svend Birkelund and Gunna Christiansen. Development of real-time PCR for detection of *Mycoplasma hominis* *BMC Microbiology* 2004, 4:35 doi:10.1186/1471-2180-4-35

Ayling, R. D., R. A. J. Nicholas, and K. E. Johansson. 1997. Application of the polymerase chain reaction for the routine identification of *Mycoplasma bovis*. *Vet. Rec.* 141:307–308.

Baczynska, A.; Svenstrup, H. E.; Fedder, J.; Birkelund, S.; Christiansen, G. Development of real-time PCR for detection of *Mycoplasma hominis*. *BMC Microbiol.* 4:35; 2004.

Ball, H. J., and D. Finlay. 1998. Diagnostic application of monoclonal antibody (MAb)-based sandwich ELISAs. *Methods Mol. Biol.* 104:127–132.

Bashiruddin, J. B., T. K. Taylor, and A. R. Gould. 1994. A PCR-based test for the specific identification of *Mycoplasma mycoides* subspecies *mycoides* SC. *J. Vet. Diagn. Investig.* 6:428–434.

Bell, C. A.; Uhl, J. R.; Hadfield, T. L.; David, J. C.; Meyer, R. E; Smith, T. E; Cockerill, E R.,III. Detection of *Bacillus anthracis* DNA by LightCycler PCR. *J. Clin. Microbiol.* 40:2897-2902; 2002.

Bi-Iske, G. Survey of mycoplasma infections in cell cultures and a comparison of detection methods. *Zentralbl. Bakteriol. Mikrobiol. Hyg.* 269:331-340; 1988.

Bruchmüller a; E. Pirkl b; R. Herrmann b; M. Stoermer c; H. Eichler a; H. Klüter a; P. Bugert . a Introduction of a validation concept for a PCR-based Mycoplasma detection assay. *Inf. Health* (2008)

Burrows J, Nitsche A, Bayly B, Walker E, Higgins G, Kok T: Detection and subtyping of Herpes simplex virus in clinical samples by LightCycler PCR, enzyme immunoassay and cell culture. *BMC Microbiol* 2002, 2:12.

Cassell GH, Waites KB: Mycoplasma infections. In *Infectious diseases of the fetus and Newborn infant* 4th edition. Edited by: Remington J, Klein J. Philadelphia: Saunders WB;1995:619-655.

Christiansen G, Jensen LT, Boesen T, Emmersen J, Ladefoged SA, Schiøtz LK, Birkelund S: Molecular biology of Mycoplasma. *Wien Klin Wochenschr* 1997, 109:557-561.

Dallo, S. F., and J. B. Baseman. 2000. Intracellular DNA replication and long-term survival of pathogenic mycoplasmas. *Microb. Pathog.* 29:301–309.

F. J. M. VAN KUPPEVELD, J. T. M. VAN DER LOGT, A. F. ANGULO, M. J. VAN ZOEST, W. G. V. QUINT, H. G. M. NIESTERS, J. M. D. GALAMA, AND W. J. G. MELCHERS: Genus- and Species-Specific Identification of Mycoplasmas by 16S rRNA Amplification APPLIED AND ENVIRONMENTAL MICROBIOLOGY, Aug. 1992, p. 2606-2615

Fraser CM, Gocayne JD, White O, Adams MD, Clayton RA, Fleischmann RD, Bult CJ, Kerlavage AR, Sutton G, Kelley JM: The minimal gene complement of *Mycoplasma genitalium*. *Science* 1995, 270:397-403.

Garcia, M., M. W. Jackwood, M. Head, S. Levisohn, and S. H. Kleven. 1996. Use of species-specific oligonucleotide probes to detect *Mycoplasma gallisepticum*, *M. synoviae*, and *M. iowae* PCR amplification products. *J. Vet. Diagn. Investig.* 8:56–63.

G Rawadi and O Dussurget Advances in PCR-based detection of mycoplasmas contaminating cell cultures. 10.1101/gr.4.4.199 *PCR Methods Appl.* 1995 4: 199-208

Harasawa, R.; Mizusawa, H.; Nozawa, K.; Nakagawa T.; Asada, K.; Kato, I. Detection and tentative identification of dominant mycoplasma species in cell cultures by restriction analysis of the 16S-23S rRNA intergenic spacer regions. *Res. Microbiol.* 144:489-493; 1993.

Hopert A, Uphoff CC, Wirth M, Hauser H, Drexler HG. Specificity and sensitivity of polymerase chain reaction (PCR) in comparison with other methods for the detection of mycoplasma contamination in cell lines. *J Immunol Methods* 1993; 164: 91–100.

Ishikawa Y, Kozakai T, Morita H, Saida K, Oka S, Masuo Y. Rapid detection of mycoplasma Contamination in cell cultures using SYBR Green-based real-time polymerase chain reaction, *In Vitro Cell Dev Biol Anim.* 2006 Mar-Apr;42(3-4):63-9.

J. Lindsay Oaks, Shannon L. Donahoe, Fred R. Rurangirwa, Bruce A. Rideout, Martin Gilbert, and Munir Z. Virani Identification of a Novel Mycoplasma Species from an Oriental White-Backed Vulture (*Gyps bengalensis*), *JOURNAL OF CLINICAL MICROBIOLOGY*, Dec. 2004, p. 5909–5912

J. Timenetsky, L.M. Santos, M. Buzinhani and E. Mettifogo Detection of multiple mycoplasma infection in cell cultures by PCR BPCraRz ildiaente Jcotiuoronn aol fo mf Multeipdliee aml ayncodp Blaiosmloag icinafle Rcteisoea ricnh c e(2006) 39: 907-914

Kong, F.; James, G.; Gordon, S.; Zelynski, A.; Gilbert, G. L. Species-specific PCR for identification of common contaminant mollicutes in cell culture. *Appl. Environ. Microbiol.* 67:3195-3200; 2001.

Laura McAuliffe, Richard J. Ellis, Roger D. Ayling, and Robin A. J. Nichola: Differentiation of *Mycoplasma* Species by 16S Ribosomal DNA PCR and Denaturing Gradient Gel Electrophoresis Fingerprinting *JOURNAL OF CLINICAL MICROBIOLOGY*, Oct. 2003, p. 4844–4847

LEONARD HAYFLICK AND ROBERT M. CHANOCK *Mycoplasma* Species of Man *BACTERIOLOGICAL REVIEWS*, June, 1965 Vol. 29, No. 2

March, J. B., J. Clark, and M. Brodrie. 2000. Characterization of strains of *Mycoplasma mycoides* subsp. *mycoides* small colony type isolated from recent outbreaks of contagious bovine pleuropneumonia in Botswana and Tanzania: evidence for a new biotype. *J. Clin. Microbiol.* 38:1419–1425.

Mygind T, Bircelund S, Christiansen G: DNA sequencing reveals limited heterogeneity in the 16S rRNA gene from the *rrnB* operon among five *Mycoplasma hominis* isolates. *Int J Syst Bacteriol* 1998, 48:1067-1071.

Nicholas, R. A. J. 1998. The veterinary significance of mycoplasmas. *Methods Mol. Biol.* 104:17–24.

Nicholas, R. A. J., and S. E. Baker. 1998. Recovery of mycoplasmas from animals. *Methods Mol. Biol.* 104:37–43.

Nicholas, R. A. J., L. A. Khan, B. Houshaymi, R. J. Miles, R. D. Ayling, H. Hotzel, and K. Sachse. 2002. Close phylogenetic and phenotypic relatedness between *Mycoplasma ovine/caprine* serogroup 11 and *Mycoplasma bovigenitalium*. *Syst. Appl. Microbiol.* 25:396–402.

Nicholas, R. A. J., F. G. Santini, K. M. Clark, N. M. A. Palmer, P. D. Santis, and J. B. Bashiruddin. 1996. A comparison of serological tests and gross lung pathology for the detection of contagious bovine pleuropneumonia in two groups of Italian cattle. *Vet. Rec.* 139:89–93.

Schalasta G, Arents A, Schmid M, Braun RW, Enders G: Fast and type-specific analysis of herpes simplex virus types 1 and 2 by rapid PCR and fluorescence melting-curve-analysis. *Infection* 2000, 28:85-91.

Tang, J.; Hu, M.; Lee, S.; Roblin, R. A polymerase chain reaction based method for detecting *Mycoplasma/Acholeplasma* contaminants in cell culture. *J. Microbiol. Methods.* 39:121-126; 2000.

Uphoff and HG, Detection of mycoplasma in leukemia–lymphoma cell lines using polymerase chain reaction CC Drexler, *GermanyLeukemia* (2002) 16, 289–293, Nature.

Uphoff, C. C.; Drexler, H. G. Comparative PCR analysis for detection of" mycoplasma infections in continuous cell lines. *In Vitro Cell. Dev. Biol.* 38A:79-85; 2002.

Uphoff, C. C.; Drexler, H. G. Detection of Mycoplasma Contaminations. In: Helgason, C. D.; Miller C. L. ed. *Methods in molecular biology*, Vol. 290: basic cell culture protocols, 3rd ed. Totowa, NJ: Hmnana Press; 2004:13-24.

Van Kuppeverd, E. J. M.; Johansson, K.-E.; Galama, J.M.D.; Kissing, J.; BiJlske, G.; van der Logt, J.T.M.; Melchers, W. J. G. Detection of mycoplasma contamination in cell cultures by a mycoplasma group-specific PCR. *Appl. Environ. Microbiol.* 60:149-152; 1994.

Van Kuppeveld FJ, Johansson KE, Galama JM, Kissing J, Bolske G, vander Logt JT, Melchers JW: Detection of mycoplasma contamination in cell cultures by a mycoplasma group-specific PCR. *Appl Environ Microbiol* 1994, 60:149-152.

Van Kuppeveld, F. J. M., J. T. M. van der Logt, A. F. Angulo, M. J. van Zoest, W. G. V.Quint, H. G. M. Niesters, M. D.Galama, and W. J. G. Melchers. 1992. Genus- and species-specific identification of mycoplasmas by 16S rRNA amplification. *Appl. Environ. Microbiol.* 58:2606-2615. (Author's correction, 59:655, 1993.)

Victoria J. Chalker³ and Joe Brownlie Taxonomy of the canine Mollicutes by 16S rRNA gene and 16S/23S rRNA intergenic spacer region sequence comparison, *International Journal of Systematic and Evolutionary Microbiology* (2004), 54, 537–542

Weisburg, W. G., J. G. Tully, D. L. Rose, J. P. Petzel, H. Oyaizu, D. Yang, L. Mandelco, J. Sechrest, T. G. Lawrence, J. Van Etten, et al. 1989. A phylogenetic analysis of the mycoplasmas: basis for their classification. *J. Bacteriol.* 171:6455–6467.

Wilson, W. J., C. L. Strout, T. Z. DeSantis, J. L. Stilwell, A. V. Carrano, and G. L. Andersen. 2002. Sequence-specific identification of 18 pathogenic microorganisms using microarray technology. *Mol. Cell. Probes* 16:119–127.

Wirth M, Berthold E, Grashoff M, Pfutzner H, Schubert U, Hauser H. Detection of mycoplasma Contaminations by the polymerase chain reaction. *Cytotechnology* 1994; 16: 67–77

Witkin SS, Kligman I, Grifo JA, Rosenwaks Z: Ureaplasma urealyticum and Mycoplasma hominis detected by the polymerase chain reaction in the cervixes of women undergoing in vitro fertilization: prevalence and consequences. *J Assist Reprod Genet* 1995, 12:610614.

Wittwer CT, Herrmann MG, Moss AA, Rasmussen RP: Continuous fluorescence monitoring of rapid cycle DNA amplification. *Biotechniques* 1997, 22:130-138.

YOKO ISHIKAWA, TAKAHARU KOZAKAI, HATSUE MORITA, KANAME SAIDA, SYUICHI OKA, AND YOSHINOR MASUO Rapid detection of Mycoplasma contamination in cell cultures using sybr green-based real-time polymerase chain reaction. *In Vitro Cell. Dev. Biol.--Animal* 42:63-69, March and April 2006