



ZmTech Scientific: Innovation and Development of Laboratory Accessories and Bio-Technologies

® • **ZmTech Protease Inhibitors Kit (Cat. I208051)**

The cocktail (250ul 40x) supplied in DMSO, reconstituted with 10 ml Lysis buffer to obtain a 1X working solution.

Contain: Stored at -20°C

AEBSF (Serine): 40mM + Aprotinin (Serine): 40ug/ml + Bestatin (aminopeptidase): 2mM + E64 (Cysteine): 4mM + Leupeptin (Cysteine and Trypsin-like): 4mM + Pepstatin A (Aspartic): 40uM

Solution for inhibiting the wide range of protease of aspartic, cysteine, serine, aminopeptidases and metalloproteases etc, available with and without EDTA (100mM)

® • **ZmTech Phosphatase Inhibitors Kit (Cat. I208052)** Stored at -20°C

The cocktail (250ul 40x) supplied in DMSO, reconstituted with 10 ml Lysis buffer to obtain a 1X working solution.

Contain:

AEBSF (Serine): 10mM + Aprotinin (Serine): 40ug/ml + E64 (Cysteine): 10mM + Leupeptin (Cysteine and Trypsin-like): 400uM + Pepstatin A (Aspartic): 40uM + Sodium fluoride (Acid): 40mM + Sodium Orthovanadate (alkaline and tyrosine: PTPs): 80mM + sodium pyrophosphate (Ser/Thr): 10mM

Solution for inhibiting Phosphatase of both serine/threonine (Ser/Thr) and protein tyrosine phosphatases (PTPs) as well as inhibiting the broad protease of aspartic, cysteine, serine

Contact Us,

Phone: 514-702 7702 (Sam Lee) Fax: 514-254 5356 (Montreal)

Email: ZmTech.services@hotmail.com or Sam.lee@zmtechsci.com (For ordering @ Technical supports)

Web: www.zmtechsci.com

Address: 2400-8 Aird Avenue, Place Hochelaga, Montreal, QC, Canada, H1V 2W6

ZmTech Scientific Inc. endeavors to assist clients based on the highest level of customer service, competitive pricing and customer satisfaction.

Our mission is: Convenience, Speed, Safety and Economy.