



**ZmTech Scientific**: Innovation and Development of Laboratory Accessories and Bio-Technologies

® • **ZmTech Mycoplasma Disinfection Spray (Cat.:M208045)**

**Description:**

Solution for safely disinfect all laboratory surfaces and apparatuses: hoods, incubators, benches, water reservoirs cell storage boxes and liquid nitrogen containers.

Effectively against Mycoplasma, fungi, bacteria, viruses: HBV/HIV, TBC.

**Contain:**

Stored at Room Temperature

0.05% Mynox®, 9% ethanol, 32.5% 1-propanol, 0.1% glutaraldehyde, 0.02% benzelthoniacetate.

**Usage:**

Spray directly on the surfaces and wipe to dryness with a paper towel. Allow the surface to dry prior to working.

Disposable gloves and safety glasses should be worn while handling. Read the safety data sheet carefully.

**Contact Us,**

Phone: 514-702 7702 (Sam Lee)      Fax: 514-254 5356 (Montreal)

Email: [ZmTech.services@hotmail.com](mailto:ZmTech.services@hotmail.com) or [Sam.lee@zmtechsci.com](mailto:Sam.lee@zmtechsci.com) (For ordering @ Technical supports)

Web: [www.zmtechsci.com](http://www.zmtechsci.com)

Address: 2400-8 Aird Avenue, Place Hochelaga, Montreal, QC, Canada, H1V 2W6

ZmTech Scientific Inc. endeavors to assist clients based on the highest level of customer service, competitive pricing and customer satisfaction.

Our mission is: Convenience, Speed, Safety and Economy.